

The future

S-ICD with Leadless Pacing System

In development is the EMPOWER™ Modular Pacing System, which includes a leadless pacemaker and the EMBLEM MRI™ S-ICD System, and is designed to be backwards-compatible with the EMBLEM™ S-ICD family. Whether patients with life-threatening arrhythmias subsequently develop a need for pacing or vice versa, this modular solution is designed to enable doctors to treat patients with the therapies they need, when they need them.

Concept device/technology. Not available for sale.



S-ICD Randomised Clinical Trial

PRAETORIAN⁷ (n = 850) is an ongoing prospective multicentre trial, in which patients are randomised in a 1:1 ratio, either to S-ICD or Transvenous ICD. The aim is to compare ICD-related adverse events between TV-ICD and S-ICD. Results are expected in 2019.

S-ICD programming recommendations

Ventricular cut-off rate recommendations are being trialled in the UNTOUCHED¹¹ study (n = 1,100). The rate of S-ICD inappropriate shocks will be compared to the historical TV-ICD rates in the MADIT-RIT¹³ study arms B and C. Results are expected 2020.

References

1. McLeod CJ, Boersma L, Okamura H, Friedman PA. The subcutaneous implantable cardioverter defibrillator: state-of-the-art review. *Eur Heart J*. 2015.
2. 2015 ESC Guidelines for the management of patients with ventricular arrhythmias and the prevention of sudden cardiac death – European Heart Journal doi: 10.1093/eurheartj/ehv316.
3. Bardy, G. H., W. M. Smith, et al. (2010). "An entirely subcutaneous implantable cardioverter-defibrillator." *N Engl J Med* 363(1): 36-44.
4. Weiss, et al. The Safety and Efficacy of a Totally Subcutaneous Implantable-Defibrillator. *CIRCULATION*. Vol. 128, no. 9. (August 2013.): 944-953.
5. Heart Rhythm - May 2012; Vol 9:5(S1-33) AB07-2.
6. Lambiase, et al. A worldwide experience with a totally subcutaneous ICD; Preliminary results of the EFFORTLESS S-ICD Registry. *European Heart Journal* Mar2014.
7. Olde Nordkamp, L.R. et al. Rationale and design of the PRAETORIAN trial: a Prospective, randomized comparison of subcutaneous and transvenous implantable cardioverter-defibrillator therapy. *Am. Heart.J.* 2012 May;163(5):753-760. e2. doi: 10.1016/j.ahj.2012.02.012.
8. M.C. Burke et al. Safety and Efficacy of the Totally Subcutaneous Implantable Defibrillator: 2-year Results from a Pooled Analysis of the IDE Study and EFFORTLESS Registry. *JACC* 2015.
9. Priori, SG. et al. *Eur Heart J*. 2015 Aug 29. doi:10.1093/eurheartj/ehv316. Epub ahead of print.
10. S-ICD® System Post Approval Study: <https://clinicaltrials.gov/ct2/show/NCT01736618>.
11. Understanding Outcomes With the EMBLEM™ S-ICD in Primary Prevention Patients With Low Ejection Fraction (UNTOUCHED): <https://clinicaltrials.gov/ct2/show/NCT02433379>.
12. Boersma L, et al. EFFORTLESS LBCT HRS 2016.
13. Moss AJ, Schuger C, Beck CA, et al. Reduction in Inappropriate Therapy and Mortality through ICD Programming. *N Engl J Med* 2012; 367: 2275 – 2283.

All cited trademarks are the property of their respective owners.
CAUTION: Law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information not intended for distribution or use in France.
CRM_403703_AA_JUN2016. Printed in the The Netherlands by De Budelse/Gosling.

S-ICD™ System Evolution

Boston Scientific
Advancing science for life™

The journey of the Subcutaneous ICD



Boston Scientific
Advancing science for life™

www.bostonscientific.eu

© 2016 Boston Scientific Corporation
or its affiliates. All rights reserved.

Device evolution

1ST GENERATION S-ICD

SQ-RX® S-ICD

In 2002 the proof of concept study began, and the 1st generation SQ-RX® S-ICD system was born. The S-ICD system swiftly gained market approvals in the EU (2009) and the USA (2012).



2ND GENERATION S-ICD

EMBLEM™ S-ICD

Launched in 2015, the **EMBLEM™ S-ICD** system was 20% thinner with 40% greater longevity compared to the 1st generation, and was enabled for LATITUDE™ NXT remote monitoring.



3RD GENERATION S-ICD

EMBLEM™ MRI S-ICD

The newest generation is the **EMBLEM™ MRI S-ICD** system. It is labelled for 1.5 Tesla full-body MRI scans and includes two new features, SMART Pass technology and AF Monitor™. SMART Pass will help ensure patients receive appropriate device therapy when needed, by enhancing the INSIGHT™ Algorithm, AF Monitor is a new detection tool designed to alert doctors of silent or new onset AF.




Why the S-ICD?

The Subcutaneous ICD system offers effective defibrillation against sudden cardiac death without transvenous (TV) leads. The S-ICD avoids risks for those patients who don't require pacing, and supports the greater recognition of the increasing long-term risks of endocardial leads, such as systemic infection, acute and chronic displacement, pneumothorax and lead fracture.

Patient population evolution

Patient prioritisation (as per McLeod et al, 2015)




Inclusion in the ESC Guidelines

Subcutaneous ICDs are now recommended in the 2015 ESC Guidelines as a Class IIa recommendation, and should be considered as an option for ICD patients who don't require pacing (brady, ATP, CRT).²



Clinical data evolution

The safety and efficacy of the S-ICD system is supported by almost 15 years of clinical data and experience, with over 3,000 patients enrolled in completed or ongoing trials.



Long-term data shows S-ICDs are safe and effective

The **EFFORTLESS¹²** multi-national registry includes over 985 patients, with up to five years of follow-up. The data further validates the efficacy with a 97.4% arrhythmia conversion rate, and revealed that few S-ICDs were removed due to a change of patient indication (0.5% removed for an ATP need, and 0.1% removed for bradycardia need). There have been no reports of systemic blood infections, cardiac injuries or endocarditis to date.

ZERO CARDIAC INJURIES, SYSTEMIC BLOOD INFECTIONS OR ENDOCARDITIS