

TRAPPER™

Exchange Device

Compatibility

~30% smaller shaft diameter than POBA*

Facilitates 6F guide catheter usage during Complex PCI procedures

Product	TRAPPER Exchange Device (2.0F, 0.67 mm)
Corsair™ Microcatheter, Turnpike™ LP Catheter, Caravel™ Catheter, Finecross™ Micro-cross Catheter	6F**
STINGRAY™ LP CTO Re-Entry System	6F**
THREADER™ Micro-dilatation Catheter OTW	6F**
EMERGE™ PTCA Dilatation Catheter OTW	6F**
CROSSBOSS™ CTO Crossing Catheter	7F**
SYNERGY™ Stent OTW	7F**
ROTABLATOR™ Rotational Atherectomy System, 1.5 mm	8F**


*EMERGE™ MR Balloon Catheter, (2.6F, 0.91 mm). BSC Data on File.

** 6F – 2.00 mm, 7F – 2.33 mm, 8F – 2.67 mm

**MASTER THE
COMPLEX™**

TRAPPER™

Exchange Device



Description	Order Number (GTIN)	Ref/Catalog Number
Trapper Exchange Device	05391526210154	H749 3933013 0

Caution: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings, and instructions for use can be found in the product labelling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for use or distribution outside the US, France, and Japan only. All trademarks are the property of their respective owners.
IC-466604-AA APR 2017

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
DINCAR2711EA