

DIREXION™ and FATHOM™ the complete radial solution

NEW: Angled Tips Shapes and Radial Lengths!

The market leading FATHOM-16 guidewire just got better. FATHOM-16 is now available in pre-shaped tips to help access the most difficult vessels and 200+ cm lengths to give the extra reach you need for radial access procedures.



FATHOM-16 Guidewire

- NEW! 200 and 215 cm lengths
- NEW! Pre-Shaped Tip
- Unique Slotted Nitinol Design for added torqueability

DIREXION Microcatheter

- The only 155 cm microcatheters – providing additional reach for the most distal radial cases.



Multiple Sizes for a Variety of Clinical Practice Situations

New Tip Shapes and Lengths FATHOM™ -16 Steerable Guidewires

UPN	Order Number	Total Length (cm)	Nitinol Tip Length (cm)	Distal Floppy Tip Length (cm)	Proximal/Distal O.D. (inches)*	Tip Shape
M001509120	50-912	180	25	10	0.016	Angled
M001509200	50-920	200	25	10	0.016	Straight
M001509210	50-921	200	25	10	0.016	Angled
M001509300	50-930	215	25	10	0.016	Straight
M001509310	50-931	215	25	10	0.016	Angled

FATHOM™ -16 Steerable Guidewires

M001509000	50-900	140	25	10	0.016	Straight
M001509010	50-901	140	35	20	0.016	Straight
M001509100	50-910	180	25	10	0.016	Straight
M001509110	50-911	180	35	20	0.016	Straight

FATHOM™ -14 Steerable Guidewires

M001508100	50-810	200	35	10	0.014	Straight
M001508110	50-811	200	35	10	0.014	Angled
M001508140	50-814	300	35	10	0.014	Straight
M001508150	50-815	300	35	10	0.014	Angled

DIREXION™ AND DIREXION HI-FLO™ Torqueable Microcatheter Family

DIREXION Microcatheter

UPN	Order Number	Usable Length (cm)	Tip Shape	RO Markers
M001195220	19-522	155	Straight	1
M001195250	19-525	155	Bern	1
M001195280	19-528	155	J	1
M001195310	19-531	155	Swan	1
M001195330	19-533	155	Straight	2
M001195350	19-535	155	Bern	2

DIREXION HI-FLO Microcatheter

M001195420	19-542	155	Straight	1
M001195450	19-545	155	Bern	1
M001195480	19-548	155	J	1
M001195510	19-551	155	Swan	1

DIREXION Microcatheter Pre-Loaded System with FATHOM™ -16 Guidewire

UPN	Order Number	DIREXION Usable Length (cm)	DIREXION Tip Shape	FATHOM-16 Guidewire Overall length (cm)
M001195620	19-562	155	Straight	180
M001195650	19-565	155	Bern	180
M001195670	19-567	130	Bern (2RO)	180

DIREXION Microcatheter HI-FLO Pre-Loaded System with FATHOM™ -16 Guidewire

M001195720	19-572	155	Straight	180
M001195750	19-575	155	Bern	180

DIREXION Microcatheter Pre-Loaded System with TRANSEND™ -14 Guidewire

UPN	Order Number	DIREXION Usable Length (cm)	DIREXION Tip Shape	TRANSEND-14 Guidewire Overall Length (cm)
M001195850	19-585	155	Bern	190

All DIREXION and DIREXION HI-FLO Torqueable Microcatheters are compatible with most chemotherapy agents, alcohol, and DMSO. All DIREXION and DIREXION HI-FLO Torqueable Microcatheters have dynamic burst pressures of 1200 PSI.

* 0.014" = 0.36 mm & 0.016" = 0.41 mm

FATHOM STEERABLE GUIDEWIRES & DIREXION TORQUEABLE MICROCATHETER
FATHOM and DIREXION are registered or unregistered trademarks of Boston Scientific Corporation or its affiliates.

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome. This material is not approved for use or distribution in France.

10-482511-AA-JUN2017 Printed in Germany by medicalvision.

**Boston
Scientific**

Advancing science for life™

www.bostonscientific.eu

© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
DINONC3562EA